GB: Instructions for use
FR: Mode d'emploi
DE: Gebrauchsanweisung
IT: Istruzioni per l'uso
ES: Instrucciones de uso

BacRac
innovation in back care

"One of the most significant factors in the outcome of a patients back problem is exercise. If the correct exercises are performed well and frequently, back problems become far easier to treat and prevent"

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GB

Using BacRac: We recommend you should initially use Bacrac 2-3 times a day up to 5-6 minutes per session. This can be repeated on a daily basis if required.

Morning use: Particularly if you wake up stiff and sore.

Evening use: Especially if you are sitting at work, or in the car for long periods.

After exercise: A short session after running, cycling, contact & racket sports to reverse joint compression.

Important please read

- 1 Follow these instructions and only use Bacrac in the way described.
- 2 Do NOT stop any prescribed medication/ treatment without advice from a suitably qualified medical professional.
- 3 Bacrac is not intended as an alternative to professional medical attention.
- 4 Acute or unremitting pain/symptoms should always be examined by a suitably qualified medical professional.
- 5 If during the use of Bacrac, you feel any worsening of symptoms/pain then STOP use immediately.

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GB: Find a warm, comfortable place to set up. Position a small pillow or rolled towel to support your head/neck.

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GB: Lie on your back with your knees bent and feet apart (Shoulder width).

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GB: Slide BacRac under your lower back with the larger rib sections level with your hip.

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GB: Settle onto the support and adjust to the most comfortable position. (Bear in mind that Bacrac needs to grip to be effective)

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GB: Once in position, simply rock from side to side using one leg at a time.

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GB: The range of movement necessary is small (10-2 on the clockface). Relax into the Bacrac and STOP if you feel any increase in pain.

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GB: Maintain a slow, relaxed rocking action and re-adjust your position as the spine lengthens.

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GB: After 5-6 minutes, either slide off the Bacrac or lift up and off to one side.

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Thank you

If you have any questions, problems or simply wish to share your experience, then please contact us. We'd love to hear from you.

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