ASOISM

EN - USER INSTRUCTIONS

Product Description:

The ankle brace is designed to provide support and stability, enhancing safety during physical activity. Made with breathable and lightweight material, the brace comfortably adapts to various anatomies. It is useful for injury prevention and expediting recovery. Thanks to the velcro closure system, you can adjust the size and compression level according to the user's specific needs.

Indications:

- Prevention of ankle injuries during sports or daily activities
- Support and stabilization in case of weak ankles
- Post-injury recovery

Instructions for Use:

Wearing the Brace:

- Fully open the velcro of the brace.
- Insert the foot into the brace, ensuring that the heel is well positioned in the designated space.
- Wrap the brace straps around the ankle, starting from the bottom up.

Brace Adjustment:

- Pull the velcro straps to adjust compression for desired comfort and support.
- Ensure that the brace is neither too tight nor too loose. Excessive compression may impede proper circulation, while too little may not provide the necessary support.

Activity Use:

- Use the brace during physical activity as directed.
- If excessive pain or discomfort is experienced, discontinue use and consult a physician.

Maintenance and Care:

- Hand wash the brace in cold water with a mild detergent.
- Do not wring the brace to dry it. Gently squeeze out excess water and let it air dry away from direct heat sources.
- Periodically check the integrity of the velcro and seams.

Warnings:

- Do not use the brace if you have severe circulatory problems without consulting a physician.
- In case of skin irritation, increased pain, or discomfort, discontinue use immediately and consult a healthcare professional.

Storage:

- Store the brace in a cool, dry place when not in use.
- Avoid prolonged exposure to direct sunlight.

Customer Service: For any questions or assistance, please contact us at the following:

- Website: www.fitpaddy.com
- Email: customercare@fitpaddy.com
- Manufacturer Address: FITPADDY Srl Via dell'Artigianato 17 - 20835 - Muggiò (MB) Italy

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