CSOISM

EN - USER INSTRUCTIONS

Product Description:

Calf supports are designed to provide compression and support during physical activity, enhancing performance and accelerating muscle recovery. Packaged in a set of two, these supports are made with breathable materials that promote circulation and help prevent fatigue and injuries. Ideal for athletes and regular exercisers.

Indications:

- Enhancing sports performance
- Prevention of fatigue and injuries during exercise
- Support and compression to accelerate muscle recovery

Instructions for Use:

Wearing the Supports:

- Ensure the supports are correctly oriented with the outer seams positioned on the sides of the calves.
- Slide the support onto the foot and pull it up to completely cover the calf.

Adjusting the Supports:

- Evenly distribute compression along the calf, ensuring there are no wrinkles or overly tight areas that may restrict circulation.
- The supports should be snug but should not cause pain or discomfort.

Use During Activity:

- Wear the supports during physical activities such as running, cycling, or team sports.
- If you experience increased pain or discomfort, reduce the compression or remove the supports.

Maintenance and Care:

- Hand wash or machine wash the supports in cold water in a delicate laundry bag.
- Air dry away from direct heat sources and sunlight.

Warnings:

- Do not use the supports on open wounds or irritated skin.
- Discontinue use and consult a doctor if persistent pain or discomfort occurs.

Storage:

- Store the supports in a dry, cool place.
- Avoid direct exposure to heat sources or intense sunlight to preserve elasticity and functionality.

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