ES01

EN - USER INSTRUCTIONS

Product Description:

The elbow brace is designed to provide stability and targeted compression, helpful in preventing injuries, soothing pain, and promoting recovery from injuries. Equipped with an inner silicone lining, this brace ensures perfect grip and optimal support. It is particularly suitable for the recovery and management of epicondylitis and epitrochleitis.

Indications:

- Prevention of elbow injuries during physical activity or repetitive movements
- Relief of pain caused by epicondylitis (tennis elbow) and epitrochleitis (golfer's elbow)
- Support and compression to promote recovery from elbow injuries

Instructions for Use:

Wearing the Brace:

- Open the Velcro and position the brace around the elbow, ensuring that the silicone lining is directed towards the affected area.
- Close the Velcro so that the brace is secure but comfortable, without restricting circulation.

Adjustment of the Brace:

- Adjust the Velcro to achieve the desired level of compression.
- The brace should provide sufficient support without causing excessive pain or discomfort.

Use During Activity:

- The brace can be worn during daily and sports activities to protect and stabilize the elbow.
- If pain increases or significant discomfort occurs, adjust the compression or consult a healthcare professional.

Maintenance and Care:

- Hand wash the brace with cold water and mild soap.
- Do not twist the brace to avoid damaging the material; allow it to air dry away from direct heat sources.

Warnings:

- Do not use the brace on open or irritated skin.
- In case of increased pain or irritation, discontinue use and consult a medical professional.

Storage:

- Store the brace in a cool, dry place when not in use.
- Avoid direct exposure to sunlight or heat to preserve the integrity of the brace.

Customer Service: For any questions or assistance, please contact us at the following:

- Website: www.fitpaddy.com
- Email: customercare@fitpaddy.com
- Manufacturer's Address: FITPADDY Srl Via dell'Artigianato 17 20835 - Muggiò (MB) Italy

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