WSOI

EN - USER INSTRUCTIONS

Product Description
The wrist support is designed to provide support and comfort during physical activity, allowing the level of compression to be adjusted according to the user's specific needs. This accessory features a velcro closure system that facilitates precise adjustment of compression, helping to prevent injuries while ensuring excellent mobility of wrist movements.
Indications
Provides wrist support during exercise or other daily activities. Helps prevent wrist injuries through adequate compression. Ideal for those seeking increased wrist stability during activities that require repetitive movements or loads.
Instructions for Use
Wearing the Support: Fully open both velcro straps of the support. Position the support on the wrist so that the padded part is in contact with the skin. Wrap the first velcro strap around the wrist and fasten it, ensuring it is comfortable and not too tight. Wrap the second velcro strap over the first, adjusting the compression until the desired level of support and stability is reached. Ensure the velcro is well secured and the support does not move during activity.
Adjusting the Support: Tighten both velcro straps to achieve the desired level of compression. The compression should be sufficient to provide support, but not so tight as to impede blood circulation. Ensure the support is well positioned to not limit wrist mobility and guarantee maximum comfort.
Use During Activity: Use the support during activities that put the integrity of the wrist at risk, such as weight lifting, intense workouts, or repetitive activities. If the support causes pain or discomfort, slightly loosen the straps or review the support's position on the wrist.
Maintenance and Care
Hand wash the support with cold water and mild soap. Do not twist the support to wring it out. Allow it to air dry away from direct heat sources and sunlight.
Warnings
Do not use the wrist support on open wounds or irritated skin. Discontinue use if signs of skin irritation appear or if pain increases.
Storage
Store the support in a cool, dry place when not in use. Avoid prolonged exposure to sunlight to maintain the elasticity and effectiveness of the velcro.
Customer Service
For any questions or assistance, please contact us at the following:
Website: www.fitpaddy.com
Email: customercare@fitpaddy.com
Manufacturer's address: FITPADDY Srl, Via dell'Artigianato 17, 20835 - Muggiò (MB), Italy

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