LANAFORM

WELLNESS

Physio Stim

ELECTROSTIMULATOR

User Manual

EN



6


FR

38













NL













74

Therapy Position Suggestions

Presentation

Thank you for choosing the LANAFORM electrostimulator with TENS and EMS technologies.

Thanks to its triple action, the **Physio Stim** will help you with its pain therapy, muscles stimulation and massage functions. The 12 modes available will give you an electrostimulation programme to suit your needs.

 Please read all the instructions before using your device, especially this basic safety information.

Table of contents

Presentation

- 01 Product introduction
- 02 Product intended use
- 03 Contraindication
- 04 Structure of the product
- 05 Product main features
- 06 Product description
- 07 Methods of use
- 08 Safety precautions
- 09 Product maintenance
- 10 Product scrap processing
- 11 Problems / Solutions
- 12 Meaning of symbols
- 13 Executive standards
- 14 Warranty

Appendix EMC Declaration

1. Product introduction

Physio Stim is a portable multifunctional device powered by 2 x AAA batteries. It has 6 TENS modes and 6 EMS modes. The device delivers electrical pulses through 4 electrodes placed on the skin, allowing users to enjoy the treatment they want: pain relief therapy, muscle development or muscle relaxation.

The electronic stimulatory module has the operating elements of an ON/OFF Key, Display screen, Mode Selection key and Intensity Modification keys.

The LCD display screen can show the selected mode and program, output intensity, and time remaining of an applied mode. The device is supplied with accessories: electrodes, cables to connect the electrodes to the device, and batteries. These accessories must be handled with care by a mentally competent adult.

The electrode pads comply with the biocompatibility standards ISO 10993-5 (Cytotoxicity) and ISO 10993-10 (Irritation and Sensitization) and are interchangeable.

2. Product intended use

To be used for temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, back of the neck, arm, leg, and foot due to strain from exercise or normal household work activities by applying current to stimulate nerve. It is intended to stimulate healthy muscles in order to improve and facilitate muscle performance.

Digital EMS/TENS devices belong to the category of electrostimulation devices, offering three primary functions suitable for combined use:

- 01 Electrical stimulation of nerve pathways (TENS).
- 02 Electrical stimulation of muscle tissue (EMS).
- 03 A massage effect triggered by electrical signals.

These devices work by mimicking the impulses naturally occurring in our bodies, transferring them to nerve and muscle fibres via electrodes applied to the skin. The electrical impulses, while harmless and typically painless, influence stimulation transmission to nerves, nerve centres, and muscle groups in the targeted area.

Regular application is typically required for electrostimulation to have an effect. While it doesn't replace

regular muscle training, it serves as a sensible supplementary training element.

TENS, or transcutaneous electrical nerve stimulation, involves the electrical stimulation of nerves through the skin. It is an effective non-pharmacological method for treating various types of pain with minimal side effects when administered correctly. TENS works by inhibiting pain transmission to nerve fibres and increasing endorphin secretion, thus reducing pain sensation. It is clinically tested, approved, and can be self-administered under proper guidance.

TENS is clinically proven and approved for treating various conditions including back pain, sore joints, neuralgia, musculoskeletal pain, pain from circulatory disorders, and chronic pain.

EMS is a well-established method used in sports medicine and rehabilitation to complement traditional muscle training. It is utilized in sports and fitness to enhance muscle performance and achieve desired aesthetic results. EMS applications include targeted muscle strengthening and muscle relaxation to ease tension and fatigue.

Digital EMS/TENS units, equipped with massage technology, can effectively relieve muscle tension and combat fatigue, offering programs simulating real massage sensations.

EN

3. Contraindication

Before using the Physio Stim device, it's crucial to consult a medical professional as it should not be seen as a substitute for medical advice or treatment, especially if you are experiencing pain or have a medical condition.

WARNING!

EN

To ensure safety and avoid potential health risks, it is strongly advised not to use the Physio Stim in the following circumstances:

- If you have implanted electrical devices like a pacemaker, or metal implants.
- If you use an insulin pump.
- If you have an abnormally high body temperature (e.g. $> 39^{\circ}\text{C}$).
- If you have epilepsy (suspicious or confirmed), heart issues, are pregnant or during menstruation .
- If you have cancer.
- If you are insensitive to heat or electrical stimulation and/or if you have impaired consciousness.
- **Avoid using in body areas with swelling, infection, or skin disorders, and refrain from use during sleep.**
- **Avoid using on wet skin after a shower or on areas with sweat.**
- **Avoid using while operating vehicles.**
- **Not suitable for individuals with metal allergies.**
- **After surgery, avoid using the device if strong muscle contractions could affect the healing process.**
- **Do not use after consuming alcohol.**
- **Do not use if you are connected to a high-frequency surgical device.**
- **Avoid using in the case of acute or chronic gastrointestinal tract diseases.**
- **Never use the device near the heart; avoid placing stimulation electrodes on the front ribcage, particularly over the pectoral muscles, as this could increase the risk of ventricular fibrillations and induce cardiac arrest.**

EN

- **Avoid** applying on the skeletal skull structure, on the eyes, around the mouth, throat, or larynx.
- **Avoid** using on areas with bleeding or during the healing stage after surgery.
- **Never** place the device on the neck or carotid artery area or genital area.

Also note that the long-term effects of chronic electrical stimulation to the human body are not known.

Before using the device, seek advice from your doctor if any of the following conditions apply to you:

- **Severe** health conditions, especially if you suspect or have been diagnosed with high blood pressure, blood clotting disorders, a tendency towards blood vessel blockages, or recurring cancerous growths.
- **Any** skin issues.
- **Unexplained** long-lasting pain anywhere in the body.
- **Diabetes.**
- **Sensory** problems that diminish pain sensation, such as metabolic disorders.
- **Undergoing** medical treatment.
- **If** you experience any issues related to the stimulation treatment.
- **If** you have ongoing skin irritation from prolonged stimulation in the same area.

Important!

Only use the EMS/ TENS device:

- **On** adults.
- **For** its intended purpose and as outlined in the provided instructions for use. Misuse can be hazardous.
- **For** external use.

- **With the original accessories provided, which can be reordered. Failure to do so will void the warranty.**

Precautions:

- **Always remove electrodes from the skin by pulling firmly to prevent injury in case of highly sensitive skin.**
- **Keep the device away from heat sources and do not use it close to shortwave or microwave devices (e.g., mobile phones) to avoid unexpected current peaks.**
- **Avoid exposing the device to direct sunlight or high temperatures.**
- **Protect the device from dust, dirt, and humidity.**
- **Do not immerse the device in water or other liquids.**
- **The device is suitable for self-treatment.**
- **For hygiene reasons, the electrodes may only be used on one person.**
- **Immediately stop using the device if it does not work properly, or if you feel unwell or experience pain.**
- **Turn off the unit before removing or relocating electrodes to prevent unintended stimulation.**
- **Do not alter electrodes (e.g., cutting them), as this could increase current intensity beyond safe levels.**
- **Avoid activities where unexpected reactions could be hazardous.**
- **Ensure no metallic objects come into contact with electrodes during stimulation, and remove jewellery or piercings in the treatment area.**
- **Keep the device away from children.**
- **Do not confuse electrode cables with those of**

other devices, and do not connect electrodes to other devices.

- **Do not use the device with other devices** transmitting electrical impulses into the body.
- **Avoid using near flammable substances, gases, or explosives.**
- **Use only non-rechargeable batteries and always use the same battery types.**
- **Use the device while sitting or lying down initially** to minimize the risk of vagal responses (feeling faint).
- **Avoid applying moisturizing lotions or ointments before treatment, as they can increase electrode wear and cause unexpected current peaks.**
- **This device is not intended for use by children or individuals with restricted physical, sensory, or mental capabilities, unless supervised by someone responsible for their safety or properly instructed.**
- **Replace adhesive electrodes immediately if their adhesive capability decreases. Use the device with new adhesive electrodes to prevent skin injuries. Replace electrodes after 20 uses at the latest.**

4. Structure of the product

EN

The Physio Stim consists primarily of a main unit, electrode patches and electrode wires. The primary structure and components are shown **Page 2**

 **1** and described at **point 6.**

5. Product main features

The design is crafted with ergonomic principles, resulting in a practical and portable device.

- **LF (Low Frequency) electric pulses create an effective composite energy field.**

- It offers multi-functionality and can target multiple areas for relief of pain.
- The device includes various functional modules, allowing for adjustable treatment times and intensity levels.
- It features an LCD display for easy monitoring and control.

6. Product description 1

Accessories and parts

Main unit: 1 pc - Electrode: 4 pcs - Instruction manual: 1 pc - Batteries: 2*AAA - Electrode wires: 2 pcs.

- 01 Mode display
- 02 Functional selection (TENS/EMS)
- 03 Right channel intensity
- 04 Right channel intensity selection
- 05 Battery status
- 06 Timer
- 07 Left channel intensity
- 08 Function module switch
- 09 Left channel intensity selection
- 10 On/ Off, set time
- 11 AAA batteries
- 12 Connecting cables for electrodes
- 13 Electrode patches

Product power requirement

- Power source : DC 3V (2 * AAA LR03 batteries), 60 mA.
- Device safety classes : class II type BF.

Description of product main technical parameters

- Impulse frequency: 1-200Hz .
- Pulse width: 40-220 μ s.
- Impulse shape: square wave.
- The product output does not have DC component.
- The individual electric charge at maximum output amplitude exceeds 7 μ C.
- Each pulse's maximum output energy is 300mJ.
- The peak output voltage in an open measurement does not exceed 500V.
- The device can withstand the effects of output end open circuit and short circuit without experiencing performance degradation.
- Output amplitude can be adjusted across 0-18 levels, with the minimum output at 0-20 Vpp and the maximum output at 50-120 Vpp.
- Therapeutic equipment timer options include 15 minutes, 30 minutes, 45 minutes, and 60 minutes with a tolerance of $\pm 10\%$.
- Safety classification : internal electric source class and BF equipment type.
- Dimensions :
 - Main unit:
length 113mm * width 50mm * height 18.8mm.
 - Electrodes : 50mm * 50mm.
- Product software version number : A/0.
- The device has a shelf life of 3 years. The electrode

pads can be used up to 20 times and have a shelf life of 2 years.

EN

Product environmental requirement

Optimal operation environmental requirements:

- Environment temperature: +5°C-+40°C.
- Humidity level: 15%-93%RH.
- Atmospheric conditions: 700hPa-1060hPa.

Storage environmental requirements:

- Temperature: -25°C + 55°C.
- Humidity level: 0-93% RH.
- Atmospheric conditions: 700hPa - 1060 hPa.
- Note: at least 30 minutes is required for the device to return to normal operating temperature from the minimum/maximum storage temperature.

7. Methods of use

First use

- 01 Install the 2*AAA LR03 batteries according to the positive and negative electrodes of the battery image on the back cover of the main unit.
- 02 Connect the main unit and electrode patches with the electrode wires.
- 03 Turn On/ Off.

Press the "On/ Off" key to turn on the device. The LCD screen will light up and you can start the program. To turn off the device, press the "On/ Off" button for 1-2 seconds.

Use 2

- 01 Output plug for channel 1: output of the electrical signal from channel 1 after connecting the cable and the electrodes.

EN

- 02 Output plug for channel 2: output of the electrical signal from channel 2 after connecting the cable and the electrodes.
- 03 Attach the electrodes to the areas you want to stimulate. Ensure that if there is any perspiration on your body, it has been cleaned before use.
- 04 Start: press the "On/ Off" key to turn on the device, the LCD displays the state by default (the mode of both channels is TENS 1, intensity is 0, timer is 30 min).
- 05 Use the "+" and "-" keys to adjust the intensity level; the device offers 18 adjustment levels. "+" increases the intensity, while "-" decreases it.
- 06 Press the "F/M" key to select the mode. There are 6 TENS modes and 6 EMS modes available, with the first mode of TENS and EMS being a combination mode. Long-pressing the "F/M" button for 1-2 seconds switches between TENS and EMS functions.
- 07 Adjust the timing by short-pressing the "On/ Off" button while the device is on. The available timing options are 15, 30, 45 and 60 minutes.

Note: The device switches itself off automatically if it is not used for one minute. When the unit is switched on again, the LCD screen displays the menu selection and the most recently used menu flashes.

Information regarding the positioning of the electrodes

It is fundamental to the intended success of electrostimulation applications that electrodes are positioned sensibly. We recommend that you consult your doctor to establish the ideal electrode positions for your intended application area.

The figures on the display are intended as an initial aid to help you position the electrodes.

EN

Suggested therapeutic position of electrodes: 3

TENS mode	Area of application	Possible electrode positions	Recommended time (Min)
1	Shoulder pain	1-3	30
2	Back pain	4-11	30
3	Abdominal pain	12-15	30
4	Arm pain	16-21	30
5	Thigh pain	22-24	30
6	Leg pain	25-28	30

EN

EMS mode	Area of application	Possible electrode positions	Recommended time (Min)
1-6	Used to relax muscle spasms, promote physical recovery and muscles development	1-28	30

The following applies to the selection of electrode positions:

Electrode distance:

The larger the distance between electrodes, the larger the stimulated tissue volume. This applies to the area and depth of the tissue volume. At the same time, however, the stimulation intensity of the tissue falls the further the electrodes are apart. As a result, greater distances between electrodes mean a larger tissue volume is stimulated, but less intensively. Consequently, you must increase the impulse intensity to boost stimulation.

The following guideline applies to selection of the electrode distance:

- Optimal distance: approx. 5- 15 cm.
- At distances below 5 cm, the unit primarily stimulates surface structures intensively.
- At distances in excess of 15 cm, large areas and deep structures are stimulated very weakly.

Relationship between electrodes and muscle fibre structures:

EN

Adapt the current flow direction to the fibre structure of the muscle according to the muscle layer you would like to treat. If you are targeting superficial muscles, position the electrodes in parallel to the fibre structure and if you are targeting deeper layers of tissue, position the electrodes across the fibre structure. You can do this by positioning electrodes as crosses (i.e. diagonally).

As part of pain relief treatment (TENS), it is advisable either to position the electrodes of a channel so that the area affected by the pain is between the electrodes or to position one electrode directly on the area affected by the pain and the other electrode at a minimum distance of 2-3 cm. You may use the electrodes of the second channel to simultaneously treat additional areas affected by pain or use them in conjunction with the electrodes of the first channel to restrict the area affected by pain (position electrodes opposite). In this case, we once again recommend positioning electrodes as crosses.

EN Tip for a muscle relaxation treatment: always use all four electrodes for optimum treatment.

Use the electrodes on skin that is clean and preferably free from hair and grease in order to prolong the life of the electrodes. If required, clean the skin with water and remove hair prior to treatment.

If an electrode should come loose during use, the impulse intensity of both channels is reduced to the lowest level. Apply the electrode again and reset the desired impulse intensity.

Battery use

- 01 The main unit operates on 2*AAA LR03 batteries, providing approximately 50 days of usage (based on 30 minutes of use per day).
- 02 When the LCD display on the device shows "", it indicates low battery power, signalling the need for a battery replacement.
- 03 If the product is not going to be used for an extended period, it is advisable to remove the battery to prevent battery leakage, which could lead to operational issues.
- 04 Only use the specified type of battery for the device.
- 05 Ensure correct polarity when inserting the battery (+/-).
- 06 Avoid using the battery in environments exceeding 45°C, as it may negatively impact battery performance and lifespan. Proper disposal methods for used batteries should adhere to environmental protection guidelines.

8. Safety precautions

EN

Warning signs and legends are provided to ensure the safe and proper use of the product, thereby preventing harm to users and others. Below are the warning signs, along with their meanings:

	Meaning
 Contraindication for use	It indicates that there is a risk of accident or serious injury resulting from improper use.
 Warning	This indicates the potential for accidents or serious injury due to incorrect usage.
 Attention	This indicates the potential for personal injury or damage to goods due to incorrect usage.

EN

- **Avoid** using the therapeutic apparatus alongside HF (high frequency) apparatus to prevent burns or damage to the devices.
- If the patient simultaneously uses the therapeutic apparatus and HF apparatus, the massage plate may cause burns on the devices, potentially damaging them. Using the apparatus within 1 meter of shortwave or microwave therapeutic devices may result in unstable output.
- **Placing** electrode pads close to the chest increases the risk of heart fibrillations.
- **Do not** alter this equipment without authorization from the manufacturer.
- The device should not be used by individuals who cannot take care of themselves, infants, or those who are not sensitive to external stimuli.
- **Connecting** a patient to high-frequency surgical equipment simultaneously may result in burns at the site of the stimulator electrode pads and possible damage to the stimulator.
- **Operating** close to shortwave or microwave therapy ME equipment may cause instability in the stimulator output.
- **Stimulation** should not be applied across or through the head, directly on the eyes, covering the mouth, on the front of the neck (especially the carotid sinus), or with electrode pads placed on the chest, upper back, or crossing over the heart.
- **Take** precautions to prevent inhalation or accidental swallowing of small parts (including 2*AAA LR03 battery).
- **Prevent** sharp parts from damaging the product.
- **Only** use accessories specified by the manufacturer.

Contraindication for use

Before using this product, pregnant women, individuals during menstruation, those with sensitive skin, heart conditions, abnormal blood pressure, malignant tumours, cerebrovascular conditions, acute illnesses, or anyone under medical treatment must consult a doctor.

01 Use of the device is not recommended by individuals with skin perception disturbances or those insensitive to heat.

02 Avoid use during bathing, sweating, or sleeping.

03 For patients with cerebral haemorrhage: avoid using during unstable phases; individuals with sequelae should use the product under the supervision of doctors.

04 Use of the device is contraindicated for those who have purulent inflammation, acute blood poisoning or continuous hyperpyrexia.

05 Use of the device is contraindicated for those who have acute cardiovascular and cerebrovascular diseases.

Please stop using it immediately and consult a doctor in case of feeling unwell or skin discomfort during use.

Precautions

Please refrain from using the product on areas nearest to the heart, head, eyes, front neck (especially the carotid artery), lower back, oral cavity, or genital area, as well as on areas with skin disease.	
01 Ensure to power off the device before moving or changing accessories during use to prevent strong stimulation. 02 Children and individuals unable to express their consent should not use the product. 03 If you feel unwell while using the product, discontinue use immediately and consult a doctor. 04 Always unplug the device when not in use. 05 Do not use the product with other medical electronic devices, such as pacemakers or electrocardiographs, as it may pose a danger. 06 Avoid using the product in areas with direct sunlight, high heat, flammables, electromagnetic radiation, or humidity. 07 Do not disassemble, repair, or modify the device to prevent accidents or failure. 08 Check the equipment before each use to prevent accidents caused by exposed wires or damage. 09 Clean the device regularly with a dry, soft cloth to maintain performance, as dust may affect its function. Before each use, please ensure to check if the electrode is loose. Failure to do so may adversely affect performance or lead to other issues.	

9. Product maintenance

- 01** If the device displays an invalid key message after pressing it, even though it seems to be functioning correctly, turn it off for a few seconds. Then, turn it back on. If it still doesn't function properly, inspect the device for any signs of damage or malfunction. If you feel unwell while using the product, discontinue use immediately and consult a doctor.
- 02** If the device displays normal functioning but does not emit any electrical stimulations after you turned it on then check whether the electrodes are in direct contact with the skin (check for any hair, clothes, etc. between the electrodes and the skin that may prevent the device to operate properly).
- 03** Before cleaning the device, shut it down and place the electrode wires away from electric supply socket. You can clean the device with a soft dry cloth or towel.
- 04** Clean the electrodes after each use. You can use a slightly damp soft cloth and dry the device directly after. When the electrodes are particularly dirty, you can clean using a soft cloth or towel imbibed with a bit of medicinal alcohol (alcohol concentration: 75%).
- 05** The apparatus should be placed in a dry, ventilated and well insulated place.
- 06** When the device is moved, it must be handled with care, avoiding shaking.
- 07** Make sure to inspect the condition of the batteries in the device and verify the quality of the output. If any irregularities are detected, it is important to address them promptly and take the necessary steps to carry out the necessary repairs.

- 08** Non-professionals should not attempt to disassemble the main unit to avoid the risk of electric shock or damage to the device. Such actions could result in accidents, putting you at risk.

10. Product scrap processing

The therapeutic equipment, classified as a medical device, undergoes scrap processing according to local legislation if the complete machine has aged and is no longer usable.

Advice regarding the disposal of waste:

- All the packaging is composed of materials that pose no hazard for the environment and which can be disposed of at your local sorting centre to be used as secondary raw materials. The cardboard may be disposed of in a paper recycling bin. The packaging film must be taken to your local sorting and recycling centre.
- When you have finished using the device, please dispose of it in an environmentally friendly way and in accordance with the law. Prior to disposal, please remove the battery and dispose of it at a collection point so that it can be recycled.

Please note: under no circumstances should spent batteries be disposed of with household waste.

EN

11. Problems / Solutions

The device does not operate when the ON/OFF button is pressed. How to proceed:

- 01** Ensure batteries have been inserted correctly and are in contact with the terminals.
- 02** Replace batteries, if applicable.
- 03** Contact Customer Services.

Electrodes do not adhere to the body. How to proceed:

- 01 Clean the adhesive surface of the electrodes using a damp, lint-free cloth. Replace the electrodes if they still do not adhere securely.
- 02 Clean the skin prior to any application; do not use skincare lotions or oils prior to treatment. Shaving the part of the skin where you want to place the electrodes will help them last longer.

There is no noticeable stimulation. How to proceed:

- 01 Press the ON/OFF button to interrupt the program. Check the connection cables are correctly connected to the device and the electrodes. Ensure the electrodes are in firm contact with the treatment area.
- 02 Press the ON/OFF button to restart the program.
- 03 Check electrode positions and ensure that adhesive electrodes do not overlap.
- 04 Gradually increase the impulse intensity.
- 05 The batteries are almost empty. Replace the batteries.

You have an unpleasant sensation at the electrodes. How to proceed:

- 01 The electrodes are not positioned correctly. Check their positions and re-position, if necessary.
- 02 The electrodes are worn. This may cause irritated skin as even distribution of the current across the entire area is no longer guaranteed. For this reason, replace the electrodes.

Skin in the treatment area turns red.

How to proceed:

Immediately stop treatment and wait until your skin has returned to its normal condition. If the redness is under the electrode and disappears quickly, there is no risk – this is caused by the locally stimulated, increased blood flow.

However, consult your doctor before you continue treatment if the skin irritation persists and this is accompanied by an itchy sensation or inflammation. This may be caused by an allergic reaction to the adhesive surface.

12. Meaning of symbols

Symbol	Meaning
	Product batch number tag
	Product serial number tag
	Manufacturer
	Production date
	Precautions
	EU representative
	Warning against danger

Applied part of type BF

This product is only for indoor use

Prohibitions

Enforcement of the law

Keep out of the sun

Keep away from rain

Consult local authorities to ascertain the appropriate disposal method for potentially biohazardous parts and accessories.

CE label

IP22

The dustproof and waterproof rating ensures protection against solid objects larger than 12mm and prevents harmful effects even when inclined up to 15 degrees, resisting the intrusion of water drops.

	Medical device
	This symbol is employed to guide the user to consult documentation for further details regarding the system's use or description.

13. Executive standards

The product conforms to the following standards and laws:

- 01 Medical electrical equipment - Part 1 : General requirements for basic safety and essential performance (EN 60601-1 :2006 + A1 :2013).
- 02 Medical electrical equipment - General requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment (EN60601-1-11 :2015).
- 03 Medical electrical equipment - Part 2-10 : Particular requirements for the safety of nerve and muscle stimulators (EN60601-2-10:2015) .
- 04 Medical electrical equipment - Part 1-2 : General requirements for safety - Collateral standard: Electromagnetic compatibility - Requirements and tests (EN60601-1-2:2015).

EN

Warranty

Limited warranty

LANAFORM guarantees this product against any material or manufacturing defect for a period of two

years from the date of purchase, except in the circumstances described below.

EN

The LANAFORM guarantee does not cover damage caused as a result of normal wear to this product. In addition, the guarantee covering this LANAFORM product does not cover damage caused by abusive or inappropriate or incorrect use, accidents, the use of unauthorised accessories, changes made to the product or any other circumstance, of whatever sort, that is outside LANAFORM's control.

LANAFORM may not be held liable for any type of circumstantial, indirect or specific damage.

All implicit guarantees relating to the suitability of the product are limited to a period of two years from the initial date of purchase as long as a copy of the proof of purchase can be supplied.

On receipt, LANAFORM will repair or replace your appliance at its discretion and will return it to you. The guarantee is only valid through LANAFORM's Service Centre. Any attempt to maintain this product by a person other than LANAFORM's Service Centre will render this guarantee void.

EN

Appendix EMC Declaration

The equipment is designed for use within the electromagnetic environment described below. It is the responsibility of the customer or user to ensure its operation within such an environment. The Physio Stim is appropriate for use in professional health-care settings, excluding areas with sensitive equipment or significant electromagnetic disturbances, such as RF shielded rooms in magnetic resonance imaging systems, operating rooms near active AF surgical equipment, electrophysiology laboratories,

shielded rooms, or areas where short wave therapy equipment is utilised.

- Avoid using the system near strong electric fields, electromagnetic fields (e.g., MRI scan rooms), and mobile wireless communication devices. Operating the device in an unsuitable environment could result in malfunctions or damage.
- Compliance with EMC and EMI regulations cannot be ensured when using modified cables or those that do not meet the same standards under which the equipment was validated.
- The system should not be used adjacent to or supported by other equipment. It is important to adhere to the recommendations outlined in this manual.
- Avoid using accessories, internal components, and cables other than those specified by the manufacturer, as this may lead to increased emissions or reduced electromagnetic immunity, resulting in improper operation.
- Portable RF communications equipment, including peripherals such as antenna cables and external antennas, should be used no closer than 30cm to any part of the system, including manufacturer-specified cables. Failure to do so may degrade the equipment performance.
- To ensure basic safety in relation to electromagnetic disturbances throughout its expected service life, always use the system in the specified electromagnetic environment and follow the maintenance recommendations provided in this manual.

EN

The following tables provide information on compliance of the equipment according to the standard EN 60601-1- 2:2015.

Table 1 - Compliance class

Emissions test	Compliance	
RF emissions CISPR 11	Group 1t	
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Table 2 - Compliance standards

Phenomenon	Basic standard of	
Electrostatic discharge	IEC 61000-4-2	
Radiated RF EM fields 1	IEC 61000-4-3	
Proximity fields from RF wireless communication equipment	IEC 61000-4-3	
Electrical fast transient / bursts	IEC 61000-4-4	
Conducted disturbances induced by RF fields	IEC 61000-4-6	
Rated power Magnetic field frequency	IEC 61000-4-8	

	Electromagnetic environment and guidance
	The equipment uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment
	The equipment is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

	Immunity test	Level
	± 8KV contact ± 2 KV, ± 4 KV ± 8 KV, ± 15 KV air	± 8KV contact ± 2 KV, ± 4 KV ± 8 KV, ± 15 KV air
	3V/m 80 MHz-2.7 GHz 80% AM at 1 KHz	3V/m 80 MHz-2.7 GHz 80% AM at 1 KHz
	See table	See table
	± 1KV 100 KHz frequency repetition	± 1KV 100 KHz frequency repetition
	3V 0.15 MHz - 80 MHz 6 Vm in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1 KHz	3V 0.15 MHz - 80 MHz 6 Vm in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1 KHz
	30 A/m 50 Hz or 60 Hz	30 A/m 50 Hz or 60 Hz

Table 3 - Test specifications for enclosure port immunity to RF wireless communications equipment

Frequency test (MHz)	B a n d (MHz)	Service	
385	380-390	TETRA 400	
450	430-470	GMRS 460, FRS 460	
710	704-787	LTE 13.17 Band	
745			
780			
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE 5 Band	
870			
930			
1720	1700-1990	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE 1, 2. 4. 25 Band, UMTS	
1845			
1970			
2450	2400-2570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE 7 Band	
5240	5100-5800	WLAN 802.11 a/n	
5500			
5785			

	Modulation	Maximum Power (W)	Distance (m)	Immunity test level (V/m)
	Pulse modulation 18 Hz	1.8	0.3	27
	FM $\pm$ 5 KHz deviation 1 KHz sine	2	0.3	28
	Pulse modulation 217 Hz	0.2	0.3	9
	Pulse modulation 18Hz	2	0.3	28
	Pulse modulation 217 Hz	2	0.3	28
	Pulse modulation 217 Hz	2	0.3	28
	Pulse modulation 217 Hz	0.2	0.3	9

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Shenzhen Kentro Medical Electronics Co., Ltd.
Ltd (www.CE-marking.eu) Enterprise Hub, NW Business