**NATURE &
DECOUVERTES**

LAMPE DE LUMINOTHÉRAPIE ASTRE


ASTRO ENERGY DAYLIGHT

LÁMPARA DE LUMINOTERAPIA ASTRO

LICHTTHERAPIELAMP ASTRE

LÂMPADA DE LUMINOTERAPIA ASTRO

Réf. 15201640

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PLEASE CAREFULLY READ THIS MANUAL AND KEEP IT IN A SAFE PLACE

SAFETY RELATED REMARKS

 **ATTENTION !**

- Do not use this appliance in wet surroundings (e.g. in the bathroom or near a shower). Do not immerse appliance in water or any other liquid.
- This appliance is suitable in home healthcare environment.
- Before use, ensure that all packaging materials are removed and that there is no visible damage to the appliance or accessories. When in doubt, do not use the appliance and contact your dealer or the customer service address provided.
- The adaptor contains a transformer. Do not cut off the adaptor to replace it with another plug, as this can cause a hazardous situation.
- This appliance is only intended for use on the human body.
- Close supervision is necessary when this appliance is used by, or near children or invalids.
- Children should be supervised to ensure that they do not play with the appliance.
- There is no need for any special training for this appliance.
- Keep the cord away from hot objects and open flames.
- Connect the appliance only to the mains voltage listed on the type plate.
- Always remove the mains plug and let the appliance cool off before you touch it.
- Do not touch the appliance with wet hands while it is plugged in; do not allow any water to be sprayed on to it. The appliance must be operated only when it is completely dry.
- Protect the appliance from strong impacts.
- Avoid pulling the cord when removing the plug from the socket. Instead, hold the plug.
- If the adaptor, cord or the appliance is damaged, contact customer service or the dealer, since special tools are needed to make repairs. Never try to repair the appliance yourself.
- Disconnection from the power supply network is only guaranteed when the mains part is unplugged.
- This appliance is not suitable for use in the presence of flammable anesthetic mixtures with air or with oxygen or nitrous oxide.
- This appliance requires no calibration, preventive inspection and maintenance.
- This appliance is not repairable and contains no user serviceable parts.
- The operator shall not touch the appliance and the patient simultaneously.

- The appliance is protected against unauthorized use. Do not modify this appliance without authorization of the manufacturer.
- If this appliance is modified, appropriate inspection and testing must be conducted to ensure continued safe use of the equipment.

GENERAL NOTES

IMPORTANT

- This appliance is not intended for commercial or clinical use, but only for individual, private household use.
- If you are taking medication such as pain relievers, antihypertensives or antidepressants, consult your doctor before using this appliance.
- People with retinal diseases as well as diabetics need to be examined by an ophthalmologist before beginning to use this appliance.
- Please do not use if you have an eye disease such as cataracts, glaucoma, diseases of the optic nerve in general and inflammation of the vitreous humour.
- If you have health concerns of any kind, contact your general practitioner!

EMC INFORMATION

- **WARNING:** use of this appliance adjacent to or stacked with other device should be avoided because it could result in improper operation. If such use is necessary, this appliance and the other device should be observed to verify that they are operating normally. If the appliance does not light on or flickering, it means that the appliance is not operating normally.
- **WARNING:** don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.
- **WARNING:** do not change any cable without authorization of the manufacturer. Use of cables other than those specified or provided by the manufacturer of this appliance could result in increased electromagnetic emissions or decreased electromagnetic immunity of this appliance and result in improper operation.
- **WARNING:** portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm to any part of this appliance, including power cable. Otherwise, degradation of the performance of this equipment could result.
- If the appliance does not light on or flickering due to the electromagnetic disturbances, do not use it and try to place the appliance in other location.
- **WARNING:** for maintaining the basic safety and essential performance with regard to electromagnetic disturbances for the expected service life, please precaution to take if the use location is near (e.g. less than 1,5 km from) from AM, FM or TV broadcast antennas.

REPAIRS

IMPORTANT

- Avoid exposing the screen to direct sunlight, spotlight or any other heat source.
- Do not use alcohol or ammonia base liquid to clean the screen. The cleaning information provided in the section "Cleaning".
- Please do not attempt to repair the appliance yourself. This could result in serious injuries. Failure to observe this regulation shall void the warranty.
- If the appliance is damaged in any way or it is not working, it should not be operated. Please try to fix the problem according to FAQ table. If the problem still not resolved, please contact customer service or the dealer, since special tools are needed to make repairs. Never try to repair the appliance yourself. Manufacturer will provide proper guidelines and materials to assist to service personnel in parts repair.

INSTRUCTIONS FOR USE

Congratulations on the purchase of your new "Energy Daylight". The product will offer many years of reliable usage if handled correctly. Please keep your receipt, these instructions and the original packaging safe in case of future requirement.

Please read the instructions carefully to gain the most benefit from your new "Energy Daylight", although the usage guidelines in this document do not replace the advice of your qualified medical advisor.

INTENDED USE

Medical intended use of this appliance is to treat SAD (seasonal affective disorder) or winter depression. "Energy Daylight" is intended to boost people energy level, to adjust circadian rhythms and make people feel more awake. "Energy Daylight" is safe and medically proven to relief from the winter blues.

HOW DOES THE "ENERGY DAYLIGHT" WORK?

In autumn and winter, the seasons with the least sunlight because of shorter days, the lack of sunlight and a problem with certain brain chemicals stops the hypothalamus working properly. The lack of light is thought to affect the production of the hormone melatonin (that makes you sleepy).

"Energy Daylight" uses bright lights to simulate the sunlight during the darker autumn and winter, whilst eliminating the dangerous UV rays contained in sunlight. The additional light helps your brain to reduce to production of hormone melatonin, adjust the levels of melatonin that can help the symptoms of SAD or winter depression.

SEASONAL AFFECTIVE DISORDER

Seasonal Affective Disorder (SAD) as known as winter depression, winter blues. It may cause people:

- Feel low, listless and fatigued
- Need for more sleep
- Have trouble concentrating
- Less active than normal
- Lacking in energy

Symptoms occur during the same season each year. Most people with SAD start to have symptom in winter and feel better in summer.

CONTRAINDICATION

The following conditions may be contraindication for use of this appliance:

- suffer(ed) from severe depression
- recent eye surgery or a diagnosed eye condition which your doctor advised you to avoid bright light
- on specific medication (i.e. certain antidepressants, psychotropic drugs or malaria tablets);
- have high blood pressure with vitreous hemorrhage have no benefit;
- skin sensitive to light, such as systemic lupus erythematosus

Always consult your doctor before your start using this appliance if you have any of the above conditions.

POSSIBLE SIDE EFFECTS

The following side effects you may experience:

- headaches
- eyestrain
- sleeping problems
- nausea

Measures to avoid, minimize or alleviate these side effects:

- Increase the treatment distance
- Avoid using the appliance during the late evening or before sleep
- Stop using the appliance for a few days to make the side effects disappear and then try again.

SAFETY AND PERFORMANCE STANDARDS

The "Energy Daylight" has been designed and manufactured to conform to the highest standards of safety and performance including electromagnetic compatibility (EMC). The "Energy Daylight" conforms to the applicable requirements of the following:

 This appliance meets the requirement of Medical Device Directive 93/42/EEC for medical appliances.

0123

The appliance is complied with the following emission and immunity standard, and test level:

Immunity Test standard	IEC 60601-1-2 Test Level	Compliance Level
Electrostatic discharge (ESD) IEC 61000-4-2	± 8kV contact ± 2kV, ± 4kV, ± 8kV, ± 15kV air	± 8kV contact ± 2kV, ± 4kV, ± 8kV, ± 15kV air
Electrical fast transient/burst IEC 61000-4-4	± 2kV, for power supply lines ± 1kV Signal input/ output 100 kHz repetition frequency	± 2kV, for power supply lines Not applicable 100 kHz repetition frequency
Surge IEC 61000-4-5	± 0.5kV, ± 1kV differential lines ± 0.5kV, ± 1kV, ± kV 2 common mode	± 0,5 kV, ± 1 kV differential lines Not Applicable
Voltage dips, interruptions and variations on power supply input lines IEC 61000-4-11	0 % UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle	0 % UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m - 50Hz/60Hz	30 A/m - 50Hz/60Hz
Conducted RF IEC 61000-4-6	3 V - 0,15 MHz - 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz - 80 % AM at 1 kHz	3 V - 0,15 MHz - 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz - 80 % AM at 1 kHz
Radiated RF IEC 61000-4-3	10 V/m - 80 MHz - 2,7 GHz - 80 % AM at 1 kHz	10 V/m - 80 MHz - 2,7 GHz - 80 % AM at 1 kHz

Note: UT is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic Immunity

	Test Frequency (MHz)	Band (MHz)	Service	Modulation	Modulation (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
Radiated RF IEC61000-4-3 (Test specification for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment	385	380 - 390	TETRA 400	Pulse modulation 18 Hz	1,8	0.3	27
	450	430 - 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28
	710	704 - 787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0.3	9
	745						
	780						
	810						
	870	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 217 Hz	2	0.3	28
	930						
	1720						
	1845	1700 - 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 18 Hz	2	0.3	28
	1970						
	2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0.3	28
	5240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0.3	9
5500							
5785							


Emissions Test Standard	Compliance
RF emissions CISPR 11	Group 1, Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not Applicable

SIGNS AND SYMBOLS

The following symbols appear in this manual, on the “Energy Daylight” or on its accessories:

	Read operating instructions before use		Caution: Consult accompanying documents
	OFF/ ON		Ingress Protection Rating
	Class II Equipment		Date of Manufacture
	Specifies serial number of the “Energy Daylight”		Manufacturer
	Authorized Representative in the European Community		Sensitive to humidity

DEVICE AND CONTROLS

- 1 ON/OFF Switch
- 2 Screen

- 3 Stand
- 4 AC Power Cord

OPERATION

Before use, please check the appliance for damage or faults. If you notice any damage or faults on the appliance, do not use it and contact customer services or your supplier.

1. Fold the stand out. Place the appliance on a solid surface such as table. Place the appliance at an angle to your body to avoid looking at the appliance directly.
2. Place the appliance approximately 12-45 cm from your eyes. You can read, eat, work at computer, watch TV or exercise while you use the appliance.
3. Plug in the appliance and switch on the light with higher intensity by pressing the ON/OFF button. To choose a lower intensity, press the ON/OFF button again.

Recommended times to use Energy Daylight:

Duration of light therapy	Distance (High Intensity)	Distance (Low Intensity)
Approx. 30 min	18 cm	12 cm
Approx. 1 hr	30 cm	20 cm
Approx. 2 hr	45 cm	35 cm

In effect, the closer you are to the source of light, the less time you need to use it.

Note

- Push the plug completely into the mains socket.
- Ensure that there is a socket near the place where the appliance will stand.
- Lay the mains cable in such a way that no one can trip over it.
- Do not place the appliance on the position that is difficult to disconnect it.
- Do not look directly into the light continuously! Some eye pain or mild headaches may be experienced during initial sessions but these should pass. If in doubt just lower your treatment times until you feel comfortable again.
- The appliance automatically turns OFF after 4 hours of continuous use for safety protection.

4. During the darkest period of the year, repeat use the appliance every day at the right time to achieve the result. It is usually sufficient to use 30 minutes a day, it depends on light intensity and distance. It is safe to use to use it longer, if you want it. Do not use the appliance late within 3 hours of going to bed; it could disrupt your sleep.
5. Press and hold the ON/OFF button for 2 times to turn off the appliance at any time.

CLEANING

- Before you clean the appliance, always switch it off, unplug it and let it cool off.
- Clean using a damp cloth, to which you can apply a little detergent if necessary.
- Do not touch the appliance with wet hands.
- Do not let water run into the appliance or spill water onto the appliance.
- Cleaning can be done once per month.

STORAGE

When this appliance is not being used, it should be stored in a dry location out of reach of children. Do not keep the appliance connected to a power source. Follow the storage instructions provided in the “Technical Specification” Chapter.

1. Ambient temperature between -10°C and 60°C.
2. Maximum relative humidity 90%.
3. A Atmospheric pressure between 700 and 1060 hPa.

DISPOSAL

Please note that this appliance is marked with the WEEE symbol on the left. This means that this appliance must not be disposed of together with ordinary household waste, as electrical and electronic waste must be disposed of separately.

If you have any questions, please contact your local authorities responsible for waste disposal.

FAQ

Problem	Possible causes	Solutions
Device does not light up	On/ Off button switched off	Switch on the on/off button
	No Power	Connect the main plug correctly
	No Power	The main cable is defective. Contact customer services or your retailer
	LEDs have reached the end of its service life	If you need to have the appliance repaired, contact customer service or an authorized dealer
	LEDs faulty	

TECHNICAL SPECIFICATION

- **Model:** ELE018921V1
- **Dimension (WHD):** 25 x 25 x 8 cm
- **Weight:** 600g
- **Lighting elements:** LEDs
- **Light intensity:** > 10,000 lux in 20 cm
- **Mains connection:** 220-240-, 50Hz, 0.2A
- **Operating Environment:** 0 to +40 °C, and 30% to 90% relative humidity
Atmospheric pressure: 700-1060hPa
- **Storage Environment:** -10 to +60 °C, and 30% to 90% relative humidity
Atmospheric pressure: 700-1060hPa
- **Transport Environment:** -10 to +60 °C, and 30% to 90% relative humidity
Atmospheric pressure: 700-1060hPa

Product Classification

1. Class II equipment
 2. No applied part
 3. Protection against ingress of water: IP21
 4. Not category AP / APG equipment
 5. Mode of operation: continuous operation.
- **Maximum length of power cord:** 1800mm
 - **Location of emission aperture:** Screen
 - **Treatment Area:** 600cm²
 - **Risk Group of light source equipment:** Group 1 in blue light classified by IEC 60601-2-57:2011
 - **Product Service Life:** 5 years

Maximum output of optical radiation for the “Energy Daylight”:

Radiance Output	Risk group classified by IEC 60601-2-57:2011	Max. value
EUVA: Eye UV-A	Exempt Group	1.5x10 ⁻⁴ W·m ⁻²
ES: Actinic UV skin & eye	Exempt Group	1,5 x 10 ⁻⁴ W·m ⁻²
EIR: Infrared radiation hazard exposure limits for the eye	Exempt Group	6 x 10 ⁻⁴ W·m ⁻²
EH: Thermal hazard for the skin	Exempt Group	3,5 x 10 ¹ W·m ⁻²
LB: Blue light	Risk Group 1	10 x 10 ² W·m ⁻² ·sr ⁻¹
LR: Retinal thermal	Exempt Group	2 x 10 ² W·m ⁻² ·sr ⁻¹

IMPORTANT

If this appliance is damaged in any way or it is not working and you cannot find the answer to your question in FAQ. Please contact customer service or the dealer for assistance.

The “Energy Daylight” is manufactured by:
Zéro-Plus International Ltd. - Room 1004, 10/F., Join-In Hang Sing Centre,
71-75 Container Port Road, Kwai Chung, New Territories, Hong Kong
Email : sales@zeroplus.hk

European Authorised Representative:
Globalmind Consumer Electronics GmbH
An der Strusbek 50,
22926 Ahrensburg, Germany

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Zero-Plus International Ltd. - Room 1004, 10/F., Join-In Hang Sing Centre,