**NATURE &  
DECOUVERTES**

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Tension : 5V  
Courant : 1.0

# PLEASE CAREFULLY READ THIS MANUAL AND KEEP IT IN A SAFE PLACE

## WARNING

- Non-observance of the instructions may result in personal injury or material damage (electric shock, skin burns, fire). The following safety and hazard information is not only intended to protect your health and the health of other people, it should also protect the product. For this reason, pay attention to the following safety notes and include these instructions when handing over this product to others.
- The heated belt must not be used by people who are not sensitive to heat or by vulnerable persons as they are unable to respond to overheating.
- This heated belt must not be used by very young children (0-3 years old) as they are unable to respond to overheating.
- This heated belt may not be used by young children (3-8 years old) unless the controller has been set by a parent or a responsible adult and the child has been given adequate instructions on how to use the heated belt safely.
- This heated belt may be used by children over the age of 8 and by people with reduced physical, sensory or mental skills or a lack of experience or knowledge, provided that they are supervised and have been instructed on to use the heated belt safely, and are fully aware of the consequent risks of use.
- This is NOT A TOY! Children must not play with the heated belt.
- Cleaning and user maintenance must not be performed by children unless supervised.
- This heated belt is not designed for use in hospitals.
- Avoid direct and prolonged contact with the skin. Do not place directly on bare skin.
- Never use the device with open wounds, burns, hematomas, swellings, or similar conditions.
- Never use the heat device on animals
- Do not insert needles or sharp objects into the heated belt.
- Keep the device away from sharp objects or sharp edges at any time.
- Do not use when the device is folded or bunched up.
- Do not create the product by placing objects on it during storage.
- Do not wrinkle the device. Fold the product carefully for storage.
- Let the device cool off prior to folding or rolling it.
- Store the device at a dry location if it is not in use.
- Do not use if the heated belt is wet.
- This heated belt may only be used in conjunction with the provided adapter charger.
- This device uses supply voltage. Therefore, the same safety regulations apply as for any other electric device.
- Operate the device only with the stated supply voltage.
- Inspect that all components of the device are in a faultless condition prior to each start-up of the device. In the event of a defect or a malfunction, switch the device off immediately and disconnect the power plug from the power network.
- The electrical and magnetic fields emitted by this heated belt may interfere with the function of a pacemaker. However, they are still below the limits: max.

5000V/m, magnetic field strength: max 80A/m, magnetic flux density: max. 0.1 millistela. Please consult your doctor and the manufacturer of your pacemaker before using this heated belt.

- Do not pull, twist or make sharp bends in the cables.
- Please check this heated belt frequently for signs of wear and tear or damage. If you detect such signs, if the device has been used abnormally or if it is not working, return it to the supplier before turning it on again.
- If the cables of the heated belt or the charger are damaged, they must be replaced by the manufacturer,

its customer service department or similarly qualified person to avoid any risk.

- Do not try to repair the device.
- Do not use when charging.
- When this heated belt is switched on:
  - Do not place any products on it (e.g. suitcase or washing basket)
  - Heating sources, such as hot water bottles, heating wraps or similar items, must not be placed on the heated belt.

## INTENDED USE

**Important:** This heated lumbar belt is only designed to warm up the human body.

**Safety Important:** The heated belt is equipped with a SAFETY THERMAL SWITCH. The sensor prevents the heating belt from overheating thanks to an automatic switch-off when the temperature has reached the maximum temperature set (55°C).

### CHARGING THE POWER BANK

 **Important :** The power bank must be charged prior to initial use if the device has not been used for a long period.

- To do this, insert the USB cable into a computer's USB port or USB wall socket
- Plug the USB cable into the micro USB port of the power bank
- Both lights will turn from red to blue when fully charged
- The power bank should be charged at least once per year to maintain its capacity.
- While charging, do not cover the power bank with covers, cushions or similar items.
- Only charge the power bank in a dry environment at normal room temperature.

### USING THE POWER BANK

Insert the USB cable of the heated belt into the USB port of the power bank. Then insert the charged power bank into the integrated pouch in the heated belt. The heated belt is now ready to use.

#### WARNING

If the heated belt is being used over 3 hours, we recommend that you set the lowest temperature setting to avoid overheating, which may lead to burns on the skin.

## APPLYING AROMA PADS

Put few drops of the essential oil on the provided cotton pad. Insert the cotton pad to the aroma pocket.

## SWITCHING ON

Switch on the heated belt by pressing the button  and control the temperature setting: The red LED displays HIGH temperature setting (55°C). Press twice the button , the white LED displays MEDIUM temperature setting (50°C). Press a third time the button , the blue LED displays LOW temperature setting (45°C).

### **NOTE**

When the power bank has been almost completely run down, the heated belt is switched off automatically to protect the power bank. The power bank must then be charged.

## SWITCHING OFF

Press the button  repeatedly until the LED display of the button  is off.

## CARE AND MAINTENANCE

**Composition** : 100 % Polyester

- Washing not allowed
- Bleaching not allowed
- Do not dry heated belt in tumble dryers
- Do not iron the heated belt
- No chemical cleaning

Follow the cleaning instructions on the label.

### CLEANING:

1. CAUTION: Unplug the heated belt from the electrical outlet before cleaning.
2. Do not immerse the device in water.
3. Clean the device only when it is cold.
4. Clean any stains with a damp cloth.

Cleaning and maintenance should not be done by unsupervised children.

### STORAGE:

1. When the device is not in use, store it as instructed.
2. Store the heated belt in a cool, dry place.
3. Check the device regularly.
4. Do not crease the device by placing objects on it when stored.
5. In case of malfunction please return the device to the point of sale.

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